


EINMALMESSER


Inzisionsmesser / Stab Knives



Artikelnummer [Navi.]	Modell	Winkel	Breite	Box
012148	A-15F	15°	0.10 mm	VE 5 Messer
012138	A-225F	22.5°	0.10 mm	VE 5 Messer
012146	A-30F	30°	0.10 mm	VE 5 Messer
011109	A-45F	45°	0.10 mm	VE 5 Messer

EINMALMESSER



Lamelliermesser / Crescent Knives



Artikelnummer [Navi.]	Modell	Winkel	Breite	Box
014695	CPU-22AGF	bevel up, gewinkelt 45°	2.2 mm	VE 5 Messer
011119	CPU-26AGF	bevel up, gewinkelt 45°	2.6 mm	VE 5 Messer

EINMALMESSER

Bleb Messer für Bleb Revision / Bleb Knives



Artikelnummer [Navi.]	Modell	Winkel	Breite	Box
017374	BKB-10AGF	bevel up, gewinkelt 45°	1.0 mm	VE 5 Messer
019691	BKB2-10AGF	bevel up, gewinkelt 45°	1.0 mm	VE 5 Messer

HERSTELLER:

KAI Industries Co.,Ltd.

1110 Oyana, Seki City

Gifu

501-3992 Japan

Web. www.kai-europe.com

VERTRIEB:

Teleon Surgical Vertriebs GmbH

Am Borsigturm 58



13507 Berlin - Deutschland

Tel. +49 (0)30 / 43 09 55 0 | Fax. +49 (0)30 / 43 09 55 11

Email. info@teleon-surgical.com | Web. www.teleon-surgical.com

EINMALMESSER



Phakolanzen / Slit Knives



Artikelnummer [Navi.]	Modell	Winkel	Breite	Box
012152	SLU-13AGF	bevel up, gewinkelt 45°	1.3 mm	VE 5 Messer
015148	SLU-18AGF	bevel up, gewinkelt 45°	1.8 mm	VE 5 Messer
012150	SLU-20AGF	bevel up, gewinkelt 45°	2.0 mm	VE 5 Messer
014696	SLU-22AGF	bevel up, gewinkelt 45°	2.2 mm	VE 5 Messer
015146	SLU-23AGF	bevel up, gewinkelt 45°	2.3 mm	VE 5 Messer
015147	SLU-24AGF	bevel up, gewinkelt 45°	2.4 mm	VE 5 Messer
012149	SLU-25AGF	bevel up, gewinkelt 45°	2.5 mm	VE 5 Messer
012143	SLU-265AGF	bevel up, gewinkelt 45°	2.65 mm	VE 5 Messer
012142	SLU-275AGF	bevel up, gewinkelt 45°	2.75 mm	VE 5 Messer
012141	SLU-28AGF	bevel up, gewinkelt 45°	2.8 mm	VE 5 Messer
014636	SLU-285AGF	bevel up, gewinkelt 45°	2.85 mm	VE 5 Messer
012140	SLU-30AGF	bevel up, gewinkelt 45°	3.0 mm	VE 5 Messer
011114	SLU-32AGF	bevel up, gewinkelt 45°	3.2 mm	VE 5 Messer

EINMALMESSER

Erweiterungsmesser / Implant Knives



Artikelnummer [Navi.]	Modell	Winkel	Breite	Box
012147	RNU-52AGF	bevel up, gewinkelt 45°	5.2 mm	VE 5 Messer
015577	RNU-56AGF	bevel up, gewinkelt 45°	5.6 mm	VE 5 Messer

HERSTELLER:

KAI Industries Co.,Ltd.

1110 Oyana, Seki City

Gifu

501-3992 Japan

Web. www.kai-europe.com

VERTRIEB:

Teleon Surgical Vertriebs GmbH

Am Borsigturm 58


13507 Berlin - Deutschland

Tel. +49 (0)30 / 43 09 55 0 | Fax. +49 (0)30 / 43 09 55 11

Email. info@teleon-surgical.com | Web. www.teleon-surgical.com

EINMALMESSER


Vitrektomiemesser / MVR Knives



Artikelnummer [Navi.]	Modell	Gauge	Breite / Ø /Form	Box
012157	VK-19GF	MVR 19 Gauge	1.5 mm / 1.0 mm / gerade	VE 5 Messer
012155	VK-20GF	MVR 20 Gauge	1.3 mm / 0.9 mm / gerade	VE 5 Messer
012159	VK-23GF	MVR 23 Gauge	0.75 mm / 0.6 mm / gerade	VE 5 Messer
012154	VKB-20GF	MVR 20 Gauge	1.3 mm / 0.9 mm / gebogen	VE 5 Messer

EINMALMESSER

Clear Corneal Messer / Clear Corneal Knives



Artikelnummer [Navi.]	Modell	Winkel	Breite	Box
019555	CCR-18AGF	double bevel, gewinkelt 45°	1.8 mm	VE 5 Messer
014697	CCR-22AGF	double bevel, gewinkelt 45°	2.2 mm	VE 5 Messer
019248	CCR-23AGF	double bevel, gewinkelt 45°	2.3 mm	VE 5 Messer
018234	CCR-24AGF	double bevel, gewinkelt 45°	2.4 mm	VE 5 Messer
011113	CCR-25AGF	double bevel, gewinkelt 45°	2.5 mm	VE 5 Messer
011108	CCR-265AGF	double bevel, gewinkelt 45°	2.65 mm	VE 5 Messer
019534	CCR-275AGF	double bevel, gewinkelt 45°	2.75 mm	VE 5 Messer
012153	CCR-28AGF	double bevel, gewinkelt 45°	2.8 mm	VE 5 Messer
011110	CCR-30AGF	double bevel, gewinkelt 45°	3.0 mm	VE 5 Messer
011111	CCR-32AGF	double bevel, gewinkelt 45°	3.2 mm	VE 5 Messer

HERSTELLER:

KAI Industries Co.,Ltd.

1110 Oyana, Seki City

Gifu

501-3992 Japan

Web. www.kai-europe.com

VERTRIEB:

Teleon Surgical Vertriebs GmbH

Am Borsigturm 58


13507 Berlin - Deutschland

Tel. +49 (0)30 / 43 09 55 0 | Fax. +49 (0)30 / 43 09 55 11

Email. info@teleon-surgical.com | Web. www.teleon-surgical.com

EINMALMESSER

Sicherheits-Vitrektomiemesser / SAFETY MVR Knives



Artikelnummer [Navi.]	Modell	Gauge	Breite /Form	Box
015648	SK-VKB20G	MVR 20 Gauge	0.9 mm / gebogen	VE 5 Messer

EINMALMESSER

Sideport Messer / Sideport Knives



Artikelnummer [Navi.]	Modell	Form / Winkel	Breite	Box
019621	SP-T12	gewinkelt 30°	1.2 mm	VE 5 Messer
019737	SP-S15	gerade	1.0 mm	VE 5 Messer

HERSTELLER:

KAI Industries Co.,Ltd.

1110 Oyana, Seki City

Gifu

501-3992 Japan

Web. www.kai-europe.com

VERTRIEB:

Teleon Surgical Vertriebs GmbH

Am Borsigturm 58

13507 Berlin - Deutschland

Tel. +49 (0)30 / 43 09 55 0 | Fax. +49 (0)30 / 43 09 55 11

Email. info@teleon-surgical.com | Web. www.teleon-surgical.com