


KANÜLEN



Artikelnummer [Navi.]	Modell	Beschreibung	Visual	Box
019835	M2273C-02/10 Hydrodissector (Sauter)	27G x 7/8" (0.40 x 22 mm) curved, tip flattended in vertical plane		VE 10 Kanülen
019836	M2273C-02/100 Hydrodissector (Sauter)	27G x 7/8" (0.40 x 22 mm) curved, tip flattended in vertical plane		VE 100 Kanülen
019837	M2639-02/10 Anterior chamber cannula (Rycroft)	25G x 1" (0.50 x 25 mm) angled 35°, 7.5 mm from tip		VE 10 Kanülen
019839	M2804-02/10 Irrigating Lasik cannula	25G x 7/8" (0.50 x 22 mm) formed 7 mm from tip 6x 0.15 mm ports at 2.4 and 6 mm from tip		VE 10 Kanülen
019794	M3610C-02/100 Cystotome	25G x 5/8" (0.50 x 16 mm) formed, tip angled 35°		VE 100 Kanülen
019795	M3630-02/10 Cystotome	30G x 5/8" (0.30 x 16 mm) formed		VE 10 Kanülen
019840	M3909-02/10 Viscoelastic cannula	23G x 7/8" (0.64 x 22 mm) angled 45°, 9 mm from tip		VE 10 Kanülen
019841	M3926A-02/10 Viscoelastic cannula	27G x 7/8" (0.40 x 22 mm) angled 35°, 9 mm from tip		VE 10 Kanülen
019842	M8605-02/10 Silicone tipped polisher	27G x 7/8" (0.40 x 22 mm) angled		VE 10 Kanülen

HERSTELLER:

Sterimedix Ltd.

26 Thornhill Road

North Moons Moat

Redditch, Worcestershire, B98 9ND - UK

Web. www.sterimedix.com

VERTRIEB:

Teleon Surgical Vertriebs GmbH

Am Borsigturm 58

13507 Berlin - Deutschland

Tel. +49 (0)30 / 43 09 55 0 | Fax. +49 (0)30 / 43 09 55 11

Email. info@teleon-surgical.com | Web. www.teleon-surgical.com