


UNIQUE BRILLIANCE



ACUNEX[®] VARIOMAX_{toric}

THE MULTIFOCAL TORIC PREMIUM IOL
FOR BRILLIANT, NATURAL VISION

ACUNEX[®] VARIOMAX^{toric}

WITH PRESBYOPIA AND ASTIGMATISM CORRECTION


50%
of patients
do not have
astigmatism

Correction of
presbyopia is
sufficient.



The other
50%
of patients
need **more**

ACUNEX[®]
VARIOMAX^{toric}



PRESBYOPIA


Age-related loss of the close adaptability of the eye due to switching between near and far vision.

ASTIGMATISM

Impaired vision at all distances, due to irregular curvature of the cornea.

With the multifocal toric IOL ACUNEX[®] VARIOMAX^{toric} you have the opportunity to effortlessly correct presbyopia and astigmatism in the course of a cataract treatment.

One-piece multifocal-toric posterior chamber lens for optimised depth of focus with aspherical surface and blue light filter



Product	ACUNEX® VARIOMAX TORIC AN6VMT	aspheric
Type	Foldable one-piece multifocal-toric IOL for capsular bag fixation	multifocal
Optics Size	6.0 mm	blue light filter
Overall Length	12.5 mm	glistening-free
Haptic Angulation	0°	standard-toric
Optics Design	Biconvex Sector-shaped near vision segment - anterior: +3.0D Aspherical and toric surface - posterior Spherical aberration neutral	
IOL Design	C-Loop haptic	
Material	Hydrophobic glistening-free acrylate UV absorbing and additional blue light filter	
Available Diopters	SE +10.0D to +30.0D (0.5D) cyl.: T0 +0.75D T1 +1.5D T2 +2.25D T3 +3.0D	
Refractive Index	1.54 @ 35°	
A-Constant (nominal)	119.1	
Sterilisation	Gamma sterilisation	
Storage	Storage in 0.9 % saline solution	
Recommended Incision Size	≥ 2.0 mm	

MANUFACTURER:



✉ marketing@teleon-surgical.com

🌐 www.teleon-surgical.com

Teleon Surgical B.V.

Van Rensselaerweg 4 b
6956AV Spankeren, The Netherlands