


UNIQUE BRILLIANCE



ACUNEX[®] VARIO_{toric}

THE EDOF - IOL WITH
ASTIGMATISM CORRECTION

ACUNEX[®] VARIO_{toric} WITH ASTIGMATISM CORRECTION



50%
of patients do not have astigmatism

VISION
without
astigmatism



The other
50%
of patients need **more**

VISION
with
astigmatism




ACUNEX[®]
VARIO_{toric}



With the ACUNEX[®] VARIO_{toric} you have the perfect EDoF-IOL to correct astigmatism. Available with standard cylinders ranging from 0.75D to 3.00D in 0.75D steps.

One-piece toric posterior chamber lens for extended depth of focus and high contrast sensitivity with aspherical surface and blue light filter



Product	ACUNEX® VARIO TORIC AN6VT	aspheric
Type	Foldable toric one-piece EDoF-IOL for capsular bag fixation	EDoF
Optic Size	6.0 mm	blue light filter
Overall Length	12.5 mm	glistening-free
Haptic Angulation	0°	standard toric
Optic Design	Biconvex Sector-shaped near vision segment - anterior: +1.5D Aspherical and toric surface - posterior Spherical aberration neutral	
IOL Design	C-Loop haptic	
Material	Hydrophobic glistening-free acrylate UV absorbing with additional blue light filter	
Available Diopters	SE +10.0D to +30.0D (0.5D) cyl.: T0 +0.75D T1 +1.5D T2 +2.25D T3 +3.0D	
Refractive Index	1.54 @ 35°	
A-Constant (nominal)	119.1	
Sterilisation	Gamma sterilisation	
Storage	Storage in 0.9 % saline solution	
Recommended Incision Size	≥ 2.0 mm	

MANUFACTURER:



✉ marketing@teleon-surgical.com

🌐 www.teleon-surgical.com

Teleon Surgical B.V.

Van Rensselaerweg 4 b
6956AV Spankeren, The Netherlands