

LENTIS® | FEMTIS®

ACUNEX®



Manufacturer:

Viscoject Bio 2.2	Accuject 2.1-BL	Ergotouch-Fly 2.0
LP604340C	LP604595	ET1-B-GMS-20-CE
Medical AG	Medical AG	ODC Industries
Low viscosity dispersive OVD	Low viscosity dispersive OVD	Low viscosity dispersive OVD

MONOFOCAL	LENTIS® L-302-1	±0.0D to +35.0D	Yes	No	Yes
	LENTIS® L-312	-10.0D to +35.0D	Yes	No	Yes
	LENTIS® LS-312Y	±0.0D to +35.0D	Yes	No	Yes
	LENTIS® L-303	±0.0D to +35.0D	Yes	No	No
	LENTIS® L-313	-10.0D to +35.0D	Yes	No	No
	LENTIS® LS-313Y	±0.0D to +35.0D	Yes	No	No
	LENTIS® L-323	±0.0D to +35.0D	Yes	No	No
	LENTIS® LU-814VR	±0.0D to +35.0D	No	Yes	No
	FEMTIS® FB-313	+15.0D to +30.0D	Yes	No	No
	ACUNEX® AN6	+10.0D to +30.0D	No	Yes	No
Enhanced Mono.	LENTIS® QUANTUM L-333	+10.0D to +30.0D	Yes	No	No
	ACUNEX® QUANTUM AN6Q	+10.0D to +30.0D	No	Yes	No
Toric	LENTIS® QUANTUM ^{toric} LS-333 T0-T3	+10.0D to +30.0D	Yes	No	No
TORIC	LENTIS® T _{plus} LS-313 T0-T6	+10.0D to +30.0D	Yes	No	No
	LENTIS® T _{plus} ^x LU-313 T TY	-10.0D to +36.0D	Yes	No	No
	LENTIS® T _{plus} LU-814 T	-10.0D to +36.0D	No	Yes	No
EDOF	LENTIS® Comfort LS-313 MF15	-10.0D to +36.0D	Yes	No	No
	FEMTIS® Comfort FB-313 MF15	+15.0D to +30.0D	Yes	No	No
	LENTIS® Comfort LU-814 MF15	±0.0D to +36.0D	No	Yes	No
	ACUNEX® VARIO AN6V	+10.0D to +30.0D	No	Yes	No
EDOF TORIC	LENTIS® Comfort ^{toric} LS-313 MF15T0-T6	+10.0D to +30.0D	Yes	No	No
	FEMTIS® Comfort ^{toric} FB-313 MF15T0-T3	+15.0D to +30.0D	Yes	No	No
	ACUNEX® VARIOTORIC AN6VTO-T3	+10.0D to +30.0D	No	Yes	No
MULTIFOCAL	LENTIS® M _{plus} ^{toric} LS-313 MF30T0-T6	+10.0D to +30.0D	Yes	No	No
	LENTIS® M _{plus} LS-313 MF20	-10.0D to +36.0D	Yes	No	No
	LENTIS® M _{plus} LS-313 MF30	-10.0D to +36.0D	Yes	No	No
	LENTIS® M _{plus} ^x LS-313 MF30	-10.0D to +36.0D	Yes	No	No
	LENTIS® M _{plus} ^x LS-313 MFX30	-10.0D to +36.0D	Yes	No	No
	LENTIS® M _{plus} LU-814 MF20	±0.0D to +36.0D	No	Yes	No
	LENTIS® M _{plus} LU-814 MF30	±0.0D to +36.0D	No	Yes	No
	FEMTIS® M _{plus} FB-313 MF30	+15.0D to +30.0D	Yes	No	No
	ACUNEX® VARIOMAX AN6VM	+10.0D to +30.0D	No	Yes	No
MULTIFOCAL TORIC	LENTIS® M _{plus} ^{toric} LU-313 MF15T TY	±0.0D to +36.0D	Yes	No	No
	LENTIS® M _{plus} ^{toric} LU-313 MF20T TY	±0.0D to +36.0D	Yes	No	No
	LENTIS® M _{plus} ^{toric} LU-313 MF30T TY	±0.0D to +36.0D	Yes	No	No
	LENTIS® M _{plus} ^x ^{toric} LU-313 MF30T TY	±0.0D to +36.0D	Yes	No	No
	LENTIS® M _{plus} ^x ^{toric} LU-313 MFX30T TY	±0.0D to +36.0D	Yes	No	No
	FEMTIS® M _{plus} ^{toric} FB-313 MF30T0-T3	+15.0D to +30.0D	Yes	No	No
	ACUNEX® VARIOMAXTORIC AN6VMT0-T3	+10.0D to +30.0D	No	Yes	No

QF2219v14 - Revision: 09.01.2026