
TELEON®



VISCOELASTIC SOLUTIONS

HPMC (hydroxypropylmethylcellulose)

Low viscosity – more fluid (flowable) than Sodium Hyaluronate – dispersive

Product	my VISC Visco 2.0% (single-box)	my VISC Visco 2.0% (multi-box)
Concentration	2.0% hydroxypropylmethylcellulose	2.0% hydroxypropylmethylcellulose
pH-Value	6.8 - 7.5	6.8 - 7.5
Endotoxin	n.a.	n.a.
Osmolality	270 - 400 mOsmol/kg	270 - 400 mOsmol/kg
Viscosity [20°C]	4,000 - 5,000 mPa·s	4,000 - 5,000 mPa·s
Molecular Weight	70,000 - 150,000 Dalton	70,000 - 150,000 Dalton
Volume	2.0 ml	2.0 ml
Shelf Life	42 months	42 months
Storage	room temperature	room temperature
Cannula	23G	23G
Manufacturer	Albomed	Albomed
Box	single-box	multi-box (10 pcs.)
Article Number	Navision	Navision
Box + Cannula	**011033**	**011034**

MANUFACTURER:
Albomed GmbH
Hildebrandstr. 11
90592 Schwarzenbruck - Germany

Web. **www.albomed.de**

DISTRIBUTION:
Teleon Surgical Vertriebs GmbH
Am Borsigturm 58
13507 Berlin - Germany
Tel. +49 30 / 43 09 55 0 | Fax. +49 30 / 43 09 55 11
Email. info@teleon-surgical.com | Web. **www.teleon-surgical.com**

02-MA-F-0005v2