
TELEON®



# VISCOELASTIC SOLUTIONS

## Sodium Hyaluronate

Cohesive

| Product | MINIVISC plus 1.4% |
|---|---|
| **Concentration** | 1.4 % Sodium Hyaluronate |
| **pH-Value** | 6.8 - 7.6 |
| **Endotoxin** | ≤ 0.2 EU/ml |
| **Osmolality** | 200 - 400 mOsmol/kg |
| **Viscosity [20°C]** | 323,000 mPa·s |
| **Molecular Weight** | 3-4 Mio. Dalton |
| **Volume** | 1.0 ml |
| **Shelf Life** | 36 months |
| **Storage** | 2 °C - 25 °C |
| **Cannula** | 27G |
| **Manufacturer** | Bohus |
| **Box** | single-box |

| Article Number | Navision |
|---|---|
| Box + Cannula | **019832** |


**MANUFACTURER:**
**Bohus BioTech AB**
Trädgårdsgatan 4
45231 Strömstad - Sweden

Web. **www.bohusbiotech.com**

**DISTRIBUTION:**
**Teleon Surgical Vertriebs GmbH**
Am Borsigturm 58
13507 Berlin - Germany
Tel. +49 30 / 43 09 55 0 | Fax. +49 30 / 43 09 55 11
Email. info@teleon-surgical.com | Web. **www.teleon-surgical.com**


02-MA-F-0005v2