
TELEON®



VISCOELASTIC SOLUTIONS

Sodium Hyaluronate

Medium viscosity - dispersive

Product	myVISC Bio 1.0%	myVISC Bio 1.4%
Concentration	1.0 % Sodium Hyaluronate	1.4 % Sodium Hyaluronate
pH-Value	6.8 - 7.4	6.8 - 7.4
Endotoxin	n.a.	n.a.
Osmolality	270 - 400 mOsmol/kg	270 - 400 mOsmol/kg
Viscosity [20°C]	20,000 mPa·s	30,000 mPa·s
Molecular Weight	1.2 - 2.0 Mio. Dalton	1.2 - 2.0 Mio. Dalton
Volume	1.0 ml	1.0 ml
Shelf Life	42 months	42 months
Storage	2 °C - 25 °C	2 °C - 25 °C
Cannula	27G	27G
Manufacturer	Albomed	Albomed
Box	single-box	single-box
Article Number	Navision	Navision
Box + Cannula	**011035**	**011031**

MANUFACTURER:
Albomed GmbH
Hildebrandstr. 11
90592 Schwarzenbruck - Germany

Web. **www.albomed.de**

DISTRIBUTION:
Teleon Surgical Vertriebs GmbH
Am Borsigturm 58
13507 Berlin - Germany
Tel. +49 30 / 43 09 55 0 | Fax. +49 30 / 43 09 55 11
Email. info@teleon-surgical.com | Web. **www.teleon-surgical.com**

02-MA-F-0005v2